


Pneumoperitoneum Insufflation Needle



Features:

- Non-slip ergonomically designed handle.
- Sharp 14 gauge needle.
- Easy insertion.
- Blunt tip.
- Spring activated stylet with visual indicator.
- High flow on/off stopcock.
- Clear handle for easy hang drop test.
- Disposable, 10 units per box.

Catalog # IN-12, 12 CM Length

IN-15, 15 CM Length

Part # 68017-01M

68017-02M

Ordering Instructions



To place an order or to receive more information contact:
Ranfac Customer Service
800.272.6322
Monday – Friday
8:30 – 5:00 E.T.

Fax your orders to our 24 hour fax line:
508.584.8588



ISO 9001

Mail your orders to:
Ranfac® Corporation
PO Box 635
Avon, Massachusetts
02322-0635

Precision instruments for specialized applications

Ranfac® Corp. | 30 Doherty Avenue | P.O. Box 635 | Avon, Massachusetts | 02322-0635 | USA

Call toll free 1.800.2-RANFAC (1.800.272.6322) • Tel 508.588.4400 • Fax 508.584.8588 • www.ranfac.com • Email: info@ranfac.com

Since 1888

09516-01P Rev - 8/03