


● Puncture Closure Device

- Percutaneous Suturing
- Ligating Abdominal Wall Bleeders
- Trocar Closure without Extending the Incision
- Tacking up Hernia Mesh Prior to Final Placement



- 14-gauge outside diameter.
- Enables the surgeon to rapidly close each trocar site securely and safely.
- Simple and efficient one-hand operation.
- Accommodates any size suture, any size trocar.
- Cost-effective design permits multiple trocar site closure with a single device per case.
- Disposable, 10 units per box.

*Catalog #RSG-14
Part #76001-01M*

Ordering Instructions



To place an order or to receive more information contact:
Ranfac Customer Service
800.272.6322
Monday – Friday
8:30 – 5:00 E.T.

Fax your orders to our 24 hour fax line:
508.584.8588



ISO 9001

Mail your orders to:
Ranfac® Corporation
PO Box 635
Avon, Massachusetts
02322-0635

Precision instruments for specialized applications

Ranfac® Corp. | 30 Doherty Avenue | P.O. Box 635 | Avon, Massachusetts | 02322-0635 | USA

Call toll free 1.800.2-RANFAC (1.800.272.6322) • Tel 508.588.4400 • Fax 508.584.8588 • www.ranfac.com • Email: info@ranfac.com

Since 1888